

Your preferred supplier of Laboratory Heat Exchangers. Tubular Heat Exchangers, Plate Heat Exchangers or Scraped Surface Heat Exchangers.

## **Clean Bench**



Powerpoint International ++81 (0) 48 229 7830

## **Clean Bench**



## **Clean Bench**

The laminar flow Clean Bench is the most effective method of ensuring a clean zone around the critical work area. The work area can be kept at class 100 utilizing horizontal or vertical laminar flow.

·Filtered Clean Air to class 100.

·Buffer Tank ensures product never enters outside atmosphere.

•Foot switch operated filling valve. •Sampling valve inside Clean Zone.

·High Efficiency Particle Air HEPA filter and fan.









2-5-10 Hikawacho, Toda-Shi, Saitama-Ken 335-0027, Japan Tel +81 (0) 48 229 7830 Fax +81 (0) 48 229 7831 E-mail : sales@pp2k.net